CARRIERE®

*S*LX3D[™]



PERFORMANCE REVOLUTION

Designed for you, your staff & your patients with no compromises.





CARRIERE SLX3D





DO THE TWIST

With a feather touch and slight twist of the instrument, the slide is propelled to the open and locked position.



COVER IT UP

Full M/D slot coverage provides unmatched rotational control, to help ensure that you finish faster.



FEEL THIS

The smoothest brackets we've ever made.
Oh so incredibly smooth, oh so comfortable.

UNMATCHED CLINICAL CONTROL

Introducing the Carriere SLX 3D Bracket System

The dedication to practice efficiencies and to advancing the state-of-the-art of patient care will be evident when you experience the beauty, comfort, and the extraordinary performance of the all-new $Carriere^{\otimes} SLX^{\otimes} 3D^{\mathbb{T}} System$.



Don't bend wire, like a pro

With proven and precise torque values, in/out thicknesses and slot depths, the *SLX 3D* Bracket is designed to ensure less wire bending so you can spend more time doing anything but wire bending.



MEDIUM Bicuspids

NARROW Lower & upper laterals

One size doesn't fit all

The **SLX 3D Bracket System** has narrow brackets for narrow teeth, and wide brackets for wide teeth - because size matters when it comes to rapidly securing rotational control and ideal tooth positioning.



Guard rails aren't just for highways

Our exclusive Adhesive *Guard Rail* $^{\text{T}}$ **Technology** (AGR) directs excess adhesive to the sides for fast, easy clean up.



CLEAR

retention

Mushroom shaped pylons for secure



Easy on, easy off

Our compound contour bases are designed to "fit like a glove", while the beveled base edges aid in quick and clean removal.





Just do the twist

Our proprietary new *EZ Twist*® feature brings a new level of staff joy and patient comfort to your practice. For fast, comfortable wire changes, **JUST DO THE TWIST!**

A REVOLUTIONARY 3-WIRE SEQUENCE

SAGITTALFIRST

Integral to the **SAGITTAL FIRST** Treatment Philosophy is the use of the **Carriere MOTION 3D** Appliance to correct AP discrepancies and stimulate initial tooth movement.

This fosters easier transitions into larger wires earlier in treatment, less time in braces, and shorter overall treatment times.


M-SERIES" WIRES

Introducing the All-New *Carriere M-Series*[™] Wires, including .015 *Cu*[™] *Nitanium*[®] Wire and the innovative 3-Wire Sequence. *M-Series* Wires are designed to:

- Enable a simplified wire progression protocol
- Reduce the number of wires used
- Finish cases in about 12 months

3-6 MONTHS MOTION 3D APPLIANCE CLASS I PLATFORM ESTABLISHED M1 M2 M3

M-SERIES WIRE SEQUENCE - FOR METAL BRACKET SYSTEMS



M-SERIES WIRE SEQUENCE - FOR CERAMIC BRACKET SYSTEMS

Follow the 3-Wire Sequence and add the wire below prior to using the M-THREE.

WIRE NAME	DESCRIPTION	OBJECTIVE	TIME IN TREATMENT	NUMBER OF APPTS. (5-6 WEEK INTERVALS)	ITEM NUMBER
M-TWO.FIVE	.019 × .025 <i>Cu Nitanium</i> 35° with <i>EZ Stop</i> s Delivery System	For final bracket alignment prior to M-THREE wire with ceramic brackets	3-4 weeks	1 appointment	103-M2.5



EZ INVENTORY CONTROL

								SQ 38		
MAXILLARY							SLX 3D Metal		SLX 3D Clear	
	Description	Color-Code	Tooth#	Torque	Ang	Antirotation	.022 R	.022 L	.022 R	.022 L
CENTRALS	_	BLACK	U1	+17°	+5°	0°	777-UR1H-10	777-UL1H-10	776-UR1H-5	776-UL1H-5
	_	BLACK	U1	+12°	+5°	0°	777-UR1-10	777-UL1-10	776-UR1-5	776-UL1-5
LATERALS	_	PINK	U2	+12°	+8°	0°	777-UR2H-10	777-UL2H-10	776-UR2H-5	776-UL2H-5
	_	PINK	U2	+8°	+8°	0°	777-UR2-10	777-UL2-10	776-UR2-5	776-UL2-5
CUSPIDS	Hook	GREEN	U3	+7°	+9°	0°	777-UR3H-HK-10	777-UL3H-HK-10	776-UR3H-HK-5	776-UL3H-HK-5
	Hook	GREEN	U3	0°	+9°	0°	777-UR3-HK-10	777-UL3-HK-10	776-UR3-HK-5	776-UL3-HK-5
BICUSPIDS	_	PURPLE	U4 & 5's	-7°	0°	0°	777-U45-10	777-U45-10	776-U45-5	776-U45-5
	Hook	PURPLE	U4 & 5's	-7°	0°	0°	777-UR45-HK-10	777-UL45-HK-10	776-UR45-HK-5	776-UL45-HK-5
1 ST MOLARS	_	WHITE	U6's	-14°	0°	12°DO	977-UR6-HK-10	977-UL6-HK-10	976-UR6-5	976-UL6-5
2 ND MOLARS	-	GREY	U7's	-14°	0°	10°D0	977-UR7-HK-10	977-UL7-HK-10	_	_

MANDIBULAR						SLX 3D Metal		SLX 3D Clear		
	Description	Color-Code	Tooth#	Torque	Ang	Antirotation	.022 R	.022 L	.022 R	.022 L
ANTERIORS	_	YELLOW	L1&2's	-1°	0°	0°	777-L12-10	777-L12-10	776-L12-5	776-L12-5
	_	YELLOW	L1&2's	-6°	0°	0°	777-L12L-10	777-L12L-10	776-L12L-5	776-L12L-5
CUSPIDS	Hook	BLUE	L3	0°	+5°	0°	777-LR3-HK-10	777-LL3-HK-10	776-LR3-HK-5	776-LL3-HK-5
	Hook	BLUE	L3	-6°	+5°	0°	777-LR3L-HK-10	777-LL3L-HK-10	776-LR3L-HK-5	776-LL3L-HK-5
1 ST BICUSPIDS	_	GREY	L4	-12°	+2°	0°	777-LR4-10	777-LL4-10	776-LR4-5	776-LL4-5
	Hook	GREY	L4	-12°	+2°	0°	777-LR4-HK-10	777-LL4-HK-10	776-LR4-HK-5	776-LL4-HK-5
2 ND BICUSPIDS	_	RED	L5	-17°	+2°	0°	777-LR5-10	777-LL5-10	776-LR5-5	776-LL5-5
	Hook	RED	L5	-17°	+2°	0°	777-LR5-HK-10	777-LL5-HK-10	776-LR5-HK-5	776-LL5-HK-5
1 ST MOLARS	Hook	WHITE	L6's	-25°	0°	2°DO	977-LR6-HK-10	977-LL6-HK-10	_	-
2 ND MOLARS	Hook	YELLOW	L7's	-10°	0°	0°	977-LR7-HK-10	977-LL7-HK-10	_	-

CARRIERE® K	CITS (Call for availability)		SLX 3D Metal	SLX 3D Clear	
Rx	Description	QTY of sets	Total brackets	Item numbers .022	Item numbers .022
STD TORQUE	U/L 5×5 kit with hooks on 3's	1	20	777-KIT-HK	776-KIT-HK
	U/L 5×5 kit with hooks on 3's, 4's, 5's	1	20	777-KIT-HKS	776-KIT-HKS
	U/L 6×6 kit with hooks on 3's & 6's	1	24	777-KIT6-HK	776-KIT6-HK
	U/L 6×6 kit with hooks on 3's, 4's, 5's, 6's	1	24	777-KIT6-HKS	776-KIT6-HKS

SLX 3D Self-ligating Buccal Tubes incorporate all of the amazing features of **SLX 3D Brackets** to complement the revolutionary appliance system.

Another **HSO first-clear ceramic upper first molars**. Now your patients can smile all day, every day, with confidence, with the most cost-effective, complete aesthetic solution available.



THE SAGITTAL FIRST REVOLUTION

Simplify & shorten treatment with predictable results.

The **SAGITTAL FIRST** Philosophy is a ground-breaking approach that simplifies and shortens Class II and Class III treatment. It employs the **Μοτιον 3D** Appliance to treat the AP dimension at the beginning of treatment before placing brackets or aligners. By resolving the most difficult part of treatment first, you can achieve a Class I platform in 3 to 6 months, shortening total treatment time by a minimum of 6 months. You know how excited patients and parents become when you mention shorter treatment times!

PREDICTABLE

MOTION 3D Class II Appliance

INNOVATIVE

SLX 3D Bracket System

ATTACH-LESS

SLX Clear Aligners





For more information on our products and educational offerings, please contact us:

HenryScheinOrtho.com

Customer Service: **888.851.0533** | In the U.S.: **800.547.2000**Outside the U.S.: +(1) **760.448.8600** | Canada: **CERUM 800.661.9567**

To email an order: USASales@HenryScheinOrtho.com

To email an order from outside the U.S.: IntlOrders@HenryScheinOrtho.com

Ortho Organizers, Inc. 1822 Aston Ave. Carlsbad, CA 92008-7306 USA

© 2019 Ortho Organizers, Inc. All rights reserved. PN 999-345 Rev. B 08/19. US Patent 10,085.824







