





ORTHO ORGANIZERS® TRANSFORCE²™ SAGITTAL LINGUAL ARCH DEVELOPMENT PLANNER




HOW TO USE: PLACE THIS PLANNER OVER PATIENT'S STUDY MODEL. FIND THE FULLY COMPRESSED SIZE THAT BEST FITS INSIDE THE LINGUAL OUTLINE OF THE TEETH. THE PASSIVE OUTLINE SHOWS THE AMOUNT OF EXPANSION CAPABLE BY THE PREACTIVATED APPLIANCE.

NOTE: ONCE THE PLANNER HAS BEEN PRINTED, PLEASE VERIFY THE SIZE HAS NOT ALTERED. YOU MAY USE THE RULER ON THE PLANNER TO CROSS REFERENCE.



ORTHO ORGANIZERS® TRANSFORCE²™ SAGITTAL LINGUAL ARCH DEVELOPMENT PLANNER



HOW TO USE: PLACE THIS PLANNER OVER PATIENT'S STUDY MODEL. FIND THE FULLY COMPRESSED SIZE THAT BEST FITS INSIDE THE LINGUAL OUTLINE OF THE TEETH. THE PASSIVE OUTLINE SHOWS THE AMOUNT OF EXPANSION CAPABLE BY THE PREACTIVATED APPLIANCE.

NOTE: ONCE THE PLANNER HAS BEEN PRINTED, PLEASE VERIFY THE SIZE HAS NOT ALTERED. YOU MAY USE THE RULER ON THE PLANNER TO CROSS REFERENCE.

